

M-8/50, platform

DESCRIPTION

Roller platform, 8 tubes 25-30 mm diameter (50 ml tubes). Application: hybridization reactions.

SPECIFICATIONS

Capacity	8
Platform type	optional
Tube diameter	20-30 mm
Tube volume	50 ml



CAT. NUMBER

BS-010117-PK

Platform M-8/50



Multi Bio RS-24
Including PRS-26 platform
Programmable rotator

Programmable Rotator **Multi Bio RS-24** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as ...

[read more](#)