

PRS-14, platform

DESCRIPTION

14 tubes 20-30 mm diameter (50 ml tubes).

SPECIFICATIONS

Capacity	14
Platform type	optional
Tube diameter	20-30 mm
Tube volume	up to 50 ml



CAT. NUMBER

BS-010118-BK

Platform PRS-14



Multi RS-60
Including PRS-48 platform
Programmable rotator

Programmable Rotator **Multi RS-60** performs several motion types in one module. Microprocessor control allows performing not only Vertical rotation (1) of the platform, but also Reciprocal rotation (2) as well ...

[read more](#)