

PRS-22, platform

DESCRIPTION

Can accommodate 22 tubes, 10-16 mm diameter (1.5 ml - 15 ml tubes).

SPECIFICATIONS

Platform type	standard
Tube diameter	10-16 mm
Tube volume	1.5-15 ml
Capacity	22



CAT. NUMBER

BS-010117-FK

Platform PRS-22



Bio RS-24
Including PRS-22 platform
Mini-rotator

Mini-rotator **Bio RS-24** provides vertical rotation of the platform. The rotator is an ideal instrument for preventing blood coagulation in tubes and for fulfilment of procedures of biological components extraction. ...

[read more](#)