

SV-16/8, platform

DESCRIPTION

16/8/8 sockets for 1.5/0.5/0.2 ml microtest tubes (\varnothing 11/8/6 mm).

SPECIFICATIONS

| | |
|-------------|----------------|
| Diameter | 11/8/6 mm |
| Tube volume | 1.5/0.5/0.2 ml |
| Capacity | 16/8/8 |



CAT. NUMBER

BS-010210-CK

SV-16/8



MSV-3500
With all 4 platforms included
Multi Speed Vortex

Multi Speed Vortex **MSV-3500** is designed for soft or intensive mixing of reagents in different size and type plastic tubes (0.2 to 50 ml).

It is designed for operation in ...

[read more](#)