

Reagents for ELISA

DESCRIPTION

Biosan provides high-quality and professional ELISA reagents and systems. More than 400 kits available (for diagnostic and research purpose) Reagents are stable at the room temperature up to 10 days - no problems with transportation or if you forget the box overnight. All kit components are stored within the same temperature +2°- +8°C up to 12 months. Reagents are easily adapted to any open automatic ELISA processor.

All ELISA reagents are RUO certified. Biosan provide full technical and information support for your desired application installation, maintenance as well as development of full individual solution including process and instruments to meet your the most highest demands.

[Full reagent offer is available here](#)

Please contact us for more detailed information.

