

Manual extraction reagents

DESCRIPTION

Biosan provides high-quality and professional nucleic acid extraction reagents and systems. Different kit formulation ensures achieving various targets and results: liquid based, magnetic based, express method and automation-capability on popular platforms such as Qiagen Symphony and Tecan.

We offer kits for extraction from clinical samples (urine, liquor, urogenital, throat, cerfix smears, feces, blood (whole, plasma, serum)), tissues, biopsy, environmental samples, forensic samples and many more. Extracted nucleic acid are suitable for further down-stream applications such as RT-PCR, real time PCR, sequencing, amplification etc.

Products are CE and RUO certified. Biosan provide full technical and information support for your desired application installation, maintenance as well as development of full individual solution including process and instruments to meet your the most highest demands.

[Full reagent offer is available here](#)

Please contact us for more detailed information.

