

SC-24N, block

DESCRIPTION

24 x 1.5 ml microtubes

SPECIFICATIONS

Amount of 1.5 ml tubes

2



CAT. NUMBER

BS-010120-GK

SC-24N



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)