

SC-18C, block

DESCRIPTION

20 × 0.5 ml + 12 × 1.5 ml microtubes

SPECIFICATIONS

Amount of 0.5 ml tubes	20
Amount of 1.5 ml tubes	12



CAT. NUMBER

BS-010143-AK	SC-18C
--------------	--------

ACCESSORIES



TS-100C

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker TS-100C provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



TS-100C Smart

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker TS-100C Smart provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)