

# SC-24NC, block

## DESCRIPTION

24 x 1.5 ml microtubes

## SPECIFICATIONS

Amount of 1.5 ml tubes

2



## CAT. NUMBER

BS-010143-GK

SC-24NC

## ACCESSORIES



**TS-100C**

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



**TS-100C Smart**

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C Smart** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)