


B-2E, block

DESCRIPTION

Block for one deep-well plate Eppendorf® 96/1000 µl - Cat.No. 0030505204, 0030506200, 0030502205



CAT. NUMBER

BS-010159-AK	B-2E block for one deep-well plate Eppendorf® 96/1000 µl
---------------------	-------------------------------------------------------------



TS-DW
Without thermoblock
Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)