

B-2S, block

DESCRIPTION

Block for one deep-well plate Sarstedt® Megablock 96/2200 µl - Cat.No. 82.1972.002



CAT. NUMBER

BS-010159-CK

B-2S block for one deep-well plate Sarstedt®
Megablock 96/2200 µl



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)