


B-2A, block

DESCRIPTION

Block for one deep-well plate Axygen® 96/2200 µl - Cat.No. P-2ML-SQ-C



CAT. NUMBER

BS-010159-FK

B-2A Block for one deep-well plate Axygen®
96/2200 µl



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)