

R-1.5M, rotor

DESCRIPTION

Rotor for 12 x 1.5/2 ml microtest tubes

SPECIFICATIONS

Amount of 2 ml tubes

12

Amount of 1.5 ml tubes

12



CAT. NUMBER

BS-010201-AK

Rotor R-1.5M



FV-2400

Including rotors R-1.5M, R-0.5/0.2M
Minicentrifuge-Vortex Microspin

Minicentrifuge-vortex Microspin **FV-2400** is specially designed for genetic engineering research (for PCR–diagnostics experiments). Units can be used in microbiological, biochemical, clinical laboratories and industrial biotechnological laboratories.

Micro–Spin provides simultaneous mixing ...

[read more](#)