

R-0.5/0.2, rotor

DESCRIPTION

Rotor for 12 x 0.5 ml and 12 x 0.2 ml microtest tubes

SPECIFICATIONS

Amount of 0.5 ml tubes	12
Amount of 0.2 ml tubes	12



CAT. NUMBER

BS-010205-BK

Rotor R-0.5/0.2

ACCESSORIES



FVL-2400N

Including rotors R-1.5, R-0.5/0.2
Mini-centrifuge/Vortex Combi-Spin

Mini-centrifuge/vortex Combi-Spin **FVL-2400N** is specially designed for genetic engineering research (for PCR-diagnostics experiments). Units can be used in microbiological, biochemical, clinical laboratories and industrial biotechnological laboratories.

Combi-Spin provides simultaneous mixing ...

[read more](#)



MSC-3000

Including rotors R-1.5, R-0.5/0.2
Centrifuge/Vortex Multi-Spin

Centrifuge/vortex Multi-Spin **MSC-3000** is product of many years evolution of Spin-Mix-Spin technology that is intended for collecting micro volumes of reagents on the microtube's bottom (first centrifugation spin), following mixing ...

[read more](#)



MSC-6000

Including rotors R-1.5, R-0.5/0.2
Centrifuge/Vortex Multi-Spin

Centrifuge/vortex Multi-Spin **MSC-6000** is product of many years evolution of Spin-Mix-Spin technology that is intended for collecting micro volumes of reagents on the microtube's bottom (first centrifugation spin), following mixing ...

[read more](#)