

# SR-32, rotor


## DESCRIPTION

Rotor for 4 x 8-section 0,2 ml microtube strips

\* not compatible with Combi-Spins produced before 2015

## SPECIFICATIONS

Amount of 0.2 ml tubes	32
8-section strips	8



## CAT. NUMBER

BS-010205-FK

Rotor SR-32

## ACCESSORIES



#### FVL-2400N

Including rotors R-1.5, R-0.5/0.2  
Mini-centrifuge/Vortex Combi-Spin

Mini-centrifuge/vortex Combi-Spin **FVL-2400N** is specially designed for genetic engineering research (for PCR-diagnostics experiments). Units can be used in microbiological, biochemical, clinical laboratories and industrial biotechnological laboratories.

Combi-Spin provides simultaneous mixing ...

[read more](#)



#### MSC-3000

Including rotors R-1.5, R-0.5/0.2  
Centrifuge/Vortex Multi-Spin

Centrifuge/vortex Multi-Spin **MSC-3000** is product of many years evolution of Spin-Mix-Spin technology that is intended for collecting micro volumes of reagents on the microtube's bottom (first centrifugation spin), following mixing ...

[read more](#)



#### MSC-6000

Including rotors R-1.5, R-0.5/0.2  
Centrifuge/Vortex Multi-Spin

Centrifuge/vortex Multi-Spin **MSC-6000** is product of many years evolution of Spin-Mix-Spin technology that is intended for collecting micro volumes of reagents on the microtube's bottom (first centrifugation spin), following mixing ...

[read more](#)