

A-02, adapter

DESCRIPTION

Adapter 12 pieces for microtubes 0.2 ml

SPECIFICATIONS

Amount of 0.2 ml tubes

12



CAT. NUMBER

BS-010213-BK A-02 adapter - 12 pieces for microtubes 0.2 ml

ACCESSORIES



Microspin 12 Plus

Including MSL-SC lid, adapters A-02, A-05
High-speed Mini-centrifuge

High-speed mini-centrifuge **Microspin 12 Plus** is a compact desktop centrifuge designed for biomedical laboratories.

Microspin 12 Plus is used for extracting RNA/DNA samples, sedimentation of biological components, biochemical and chemical ...

[read more](#)



Microspin 12

Including MSL-SC lid, adapters A-02, A-05
High-speed Mini-centrifuge

High-speed mini-centrifuge **Microspin 12** is a compact desktop centrifuge designed for biomedical laboratories.

Microspin 12 is used for extracting RNA/DNA samples, sedimentation of biological components, biochemical and chemical analysis of ...

[read more](#)