

## RR-U, rack for rotors

### DESCRIPTION

Rack for rotors



### CAT. NUMBER

BS-010208-UK

Rack for rotors RR-U

### ACCESSORIES



**LMC-3000**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-12**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-4200R**  
Without rotor  
Laboratory Refrigerated  
Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



**LMC-56**  
Without rotor  
Laboratory Centrifuge

**LMC-56** is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



**BR-4U**  
BS-010221-AK  
Rotor

Bucket rotor without bucket adapters



**BA-2/50**  
BS-010221-CK  
Adapters

Bucket adapters for 2 tubes of 50 ml (set of 4 pcs.)



**BA-4/FT**  
BS-010221-YK  
Adapters

Bucket adapters for 4 conical bottom centrifuge tubes of 15 ml (set of 4 pcs.)



**BA-4/VT**  
BS-010221-4K  
Adapters

Bucket adapters for 4 vacutainer adapters of 2-9 ml (set of 4 pcs.). Included 8-9 ml [16 x 100] vacutainer adapters, BN-16/100-16.



**BA-4/RT**  
BS-010221-3K  
Adapters

Bucket adapters for 4 round bottom centrifuge tubes of 10 ml (set of 4 pcs.)



**BA-8/5**  
BS-010221-FK  
Adapters

Bucket adapters for 8 tubes of 5 ml (set of 4 pcs.)



**BA-14/2U**  
BS-010221-JK  
Adapters

Bucket adapters for 14 microtest tubes of 2 ml (set of 4 pcs.)



**R-6**  
BS-010208-DK  
rotor

Swing-out rotor for conical bottom centrifuge tube. Aluminium buckets.



**R-6P**  
BS-010208-XK  
Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).



**Ri-6**  
BS-010221-HK  
rotor

Swing-out rotor for conical bottom centrifuge tube. Aluminium buckets.



**Ri-6P**  
BS-010221-IK  
Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).



**RMT-24**  
BS-010221-BK  
Rotor

Fixed-angle rotor for 24 - 1.5 or 2 ml microtubes



**R-12/10**  
BS-010208-BK  
Rotor

Swing-out rotor for round bottom centrifuge tube.



**Ri-12/15**  
BS-010221-KK  
Rotor

Swing-out rotor for conical bottom centrifuge tube.



**Ri-12/10**  
BS-010221-MK  
Rotor

Swing-out rotor for round bottom centrifuge tube.



**R-24/10**  
BS-010212-JK  
Rotor

Swing-out rotor for round bottom centrifuge tube.



**Ri-24/10**  
BS-010221-LK  
Rotor

Swing-out rotor for round bottom centrifuge tube.



**R-12/15**  
BS-010208-EK  
Rotor

Swing-out rotor for conical bottom centrifuge tube.



**R-2**  
BS-010208-AK  
Rotor

Swing-out rotor for standard 96-well microtitre plates, skirted PCR plates, deepwell plates up to 45mm.



**Ri-2**  
BS-010221-GK  
Rotor

Swing-out rotor for standard 96-well microtitre plates, skirted PCR plates, deepwell plates up to 45mm.



**R-24GC**  
BS-010208-VK  
Rotor

Rotor for 24 Gel Cards for blood group serology testing (Forward Group, Reverse Group, RhD Type and 3 cell antibody screen)



**Ri-24GC**  
BS-010221-NK  
Rotor

Rotor for 24 Gel Cards for blood  
group serology testing (Forward  
Group, Reverse Group, RhD  
Type and 3 cell antibody  
screen)