

# B23-1.5, block

## DESCRIPTION

23 sockets for 1.5 ml microtest tubes depth 35 mm

## SPECIFICATIONS

Amount of 1.5 ml tubes	23
Depth	35 mm



## CAT. NUMBER

BS-010418-DK

B23-1.5



CH 3-150  
Without thermoblock  
Heating and cooling thermostat

Combitherm-2 CH 3-150 is specially designed to thermostabilise materials at temperatures from  $-3^{\circ}\text{C}$  to  $+150^{\circ}\text{C}$  according to methods of analysis. To obtain useful functionality and decrease foot-print of instruments Combitherm-2 ...

[read more](#)