


## TR-44/11, rack

### DESCRIPTION

Test tube rack for 2 / 1.5ml microtest tubes, capacity 44 tubes. WxDxH: 155x90x112mm

### SPECIFICATIONS

Capacity	44
Tube size	2 / 1.5 ml



### CAT. NUMBER

**BS-010406-JK**

Test tube rack for 2 / 1.5ml microcentrifuge tubes,  
capacity 44 tubes



**WB-4MS**  
With base BP-1 and lid  
Stirred water bath

Water bath-thermostat **WB-4MS** is designed for chemical, pharmaceutical, medical and biological laboratory research.

**WB-4MS** provides increased temperature stabilization (up to 0.1°C) due to built-in magnetic stirrer (regulated speed 250–1,000 rpm).

...

[read more](#)