

Carousel stand

DESCRIPTION

Carousel stand (rotating) for 6 pipettes

SPECIFICATIONS

Capacity	6
----------	---



CAT. NUMBER

BS-010522 Carousel stand (rotating) for 6 pipettes

ACCESSORIES



Assist
 Adjustable volume pipettes - Single channel
 pipette series

The **Assist series** pipettes are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on the model.

Multichannel ...

[read more](#)



Assist
 Adjustable volume pipettes - Single channel
 pipette series with volume lock

The **Assist series** pipettes with volume lock are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on

...

[read more](#)