

# 1-position stand

## DESCRIPTION

## SPECIFICATIONS

Capacity	1
----------	---



## CAT. NUMBER

BS-010524	1-position stand
-----------	------------------

## ACCESSORIES



**Assist**  
Adjustable volume pipettes - Single channel pipette series

The **Assist series** pipettes are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on the model.

Multichannel ...

[read more](#)



**Assist**  
Adjustable volume pipettes - Single channel pipette series with volume lock

The **Assist series** pipettes with volume lock are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on

...

[read more](#)