

## 4-position stand

### DESCRIPTION

### SPECIFICATIONS

Capacity	4
----------	---



### CAT. NUMBER

BS-010525	4-position stand
-----------	------------------

### ACCESSORIES



**Assist**  
 Adjustable volume pipettes - Single channel  
 pipette series

The **Assist series** pipettes are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on the model.

Multichannel ...

[read more](#)



**Assist**  
 Adjustable volume pipettes - Single channel  
 pipette series with volume lock

The **Assist series** pipettes with volume lock are single or 8, 12 channel variable volume pipettes designed to measure and transfer volumes from 0.1 µl to 10000 µl depending on

...

[read more](#)