

# R-24/10, Rotor

## DESCRIPTION

Swing-out rotor for round bottom centrifuge tube.

## SPECIFICATIONS

Rotor type	Swing-out
Tube dimensions (Ø × length)	16 × 105 mm
Capacity	24
Tube's volume	10–15 ml
Max. speed	4000 rpm
Max. RCF: LMC-4200R	2860 × g
Tube manufacturers	Nunc, Greiner, Greiner Bio-one, TPP etc.
Autoclavable	+




## CAT. NUMBER

BS-010212-JK

Rotor R-24/10

## ACCESSORIES



### LMC-4200R

Without rotor  
Refrigerated Laboratory  
Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



**BN-13/75-24**  
BS-010221-RK  
adapter set

For rotors R-24/10, Ri-24/10  
optional adapter set (24 pcs) -  
BN-13/75-24 for vacutainers 2-5  
ml



**BN-13/100-24**  
BS-010221-SK  
adapter set

For rotors R-24/10, Ri-24/10  
optional adapter set (24 pcs) -  
BN-13/100-24 for vacutainers 4-  
8 ml



**BN-16/100-24**  
BS-010221-TK  
adapter set

For rotors R-24/10, Ri-24/10  
optional adapter set (24 pcs) -  
BN-16/100-24 for vacutainers 8-  
10 ml



**RR-U**  
BS-010208-UK  
rack for rotors

Rack for rotors



### BN-16/90

BS-010208-SK  
bucket set for rotor R-12/10 R-  
24/10

Spare bucket set for rotor R-  
12/10, R-24/10. Plastic buckets  
(12 pcs).