

TubeSpin® Bioreactor 50 - 20

DESCRIPTION

50 ml tubes with membrane filter TPP TubeSpin® Bioreactor 50, 20 pcs.

SPECIFICATIONS

Number of tubes in set	20
Volume	50 ml



CAT. NUMBER

BS-010158-AK

50 ml tubes with membrane filter
 TubeSpin® Bioreactor 50, TPP®
 20 pcs.

ACCESSORIES



RTS-1

Including TPP TubeSpin®
Bioreactor vessels 50ml, 20pcs
Personal bioreactor

RTS-1 is personal bioreactor which utilize patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the singleuse falcon ...

[read more](#)



RTS-1C

Including TPP TubeSpin®
Bioreactor vessels 50ml, 20pcs
Personal bioreactor

RTS-1C is personal bioreactor which utilize patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the singleuse falcon ...

[read more](#)



RTS-8

Including TPP TubeSpin®
Bioreactor vessels 50ml, 20pcs
Personal multi-channel bioreactor with non-invasive real time OD measurement

RTS-8 is a personal bioreactor that utilizes patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the single-use ...

[read more](#)



RTS-8 Plus

Including TPP TubeSpin®
Bioreactor vessels 50ml, 20pcs and sterile TPP TubeSpin®
Bioreactor vessels, 50ml, with pH and O2 sensors, 10pcs
Personal multi-channel bioreactor with non-invasive real time OD pH and pO2 measurement

RTS-8 plus is a personal bioreactor that utilizes patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the ...

[read more](#)