

# F-1, shelf


## DESCRIPTION

F-1, shelf

## SPECIFICATIONS

Overall dimensions (W×D×H)

400 × 140 × 50 mm



## CAT. NUMBER

BS-040104-DK

F-1

## ACCESSORIES



### UVC/T-M-AR

With built in socket and inlet for power cords  
PCR cabinet

The Biosan **UVC/T-M-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



### UVT-B-AR

With built in socket and inlet for power cords  
PCR cabinet

The Biosan **UVT-B-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



### UVT-S-AR

BS-040107-A07  
PCR cabinet

The Biosan **UVT-S-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)