

# P-5, shelf for pipettes

## DESCRIPTION

P-5, shelf for pipettes


## SPECIFICATIONS

Overall dimensions (W×D×H)

230 × 140 × 50 mm

Capacity

5



## CAT. NUMBER

BS-040104-CK

P-5

## ACCESSORIES



### UVC/T-M-AR

With built in socket and inlet for power cords  
PCR cabinet

The Biosan **UVC/T-M-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



### UVT-B-AR

With built in socket and inlet for power cords  
PCR cabinet

The Biosan **UVT-B-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



### UVT-S-AR

BS-040107-A07  
PCR cabinet

The Biosan **UVT-S-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)