

## BN-11/30, bucket set for rotor R-6P

### DESCRIPTION

Spare bucket set for rotor R-6P. Plastic buckets (POM Kocetal) (6 pcs).

### SPECIFICATIONS

Tube holders material	POM Kocetal
Tube's volume	50 ml
Max. speed	4200 rpm
Max. temperature	150 °C
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.
Autoclavable	+



### CAT. NUMBER

BS-010208-ZK

BN-11/30

### ACCESSORIES



**LMC-3000**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-12**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-4200R**  
Without rotor  
Laboratory Refrigerated  
Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



**R-6P**  
BS-010208-XK  
Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).



**Ri-6P**  
BS-010221-IK  
Rotor

Swing-out rotor for conical bottom centrifuge tube. Plastic buckets (POM Kocetal).