

BN-17/120, bucket set for rotor R-12/15

DESCRIPTION

Spare bucket set for rotor R-12/15. Plastic buckets (12 pcs).

SPECIFICATIONS

Tube's volume	15 ml
Max. speed	4200 rpm
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.
Autoclavable	+



CAT. NUMBER

BS-010208-TK

BN-17/120

ACCESSORIES



LMC-3000
 Without rotor
 Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-12
 Without rotor
 Laboratory Centrifuge

LMC-3000 is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



LMC-4200R
 Without rotor
 Laboratory Refrigerated
 Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



Ri-12/15
 BS-010221-KK
 Rotor

Swing-out rotor for conical bottom centrifuge tube.



R-12/15
 BS-010208-EK
 Rotor

Swing-out rotor for conical bottom centrifuge tube.