

# RS2, rack for CPS-20 / CTR-6 installation

## DESCRIPTION

Rack for CPS-20 / CTR-6 installation



## CAT. NUMBER

BS-010425-HK

RS2

## ACCESSORIES



### S-Bt Smart Biotherm

Software included + RS6, rack  
with 3 shelves  
Compact CO<sub>2</sub> Incubator

S-Bt Smart Biotherm is designed for work in the areas of cell biology (operations with animal cell cultures and tissues), molecular biology (DNA/RNA reaction analysis, hybridization reactions), biotechnology (synthesis of ...

[read more](#)



### CPS-20

Without platform  
CO<sub>2</sub> Shaker

CO<sub>2</sub> Shaker CPS-20 provides regulated orbital motion of the platform and is designed for use specifically in CO<sub>2</sub> incubators. CPS-20 is specifically designed for use in harsh environments such as ...

[read more](#)



### CTR-6

BS-010174-A01  
CO<sub>2</sub> Tube roller

CO<sub>2</sub> Tube Roller CTR-6 provides regulated rocking and rolling of maximum up to 6 rollers and is designed for use specifically in CO<sub>2</sub> incubators. CTR-6 is specifically designed ...

[read more](#)