

Ri-12/10, Rotor

DESCRIPTION

Swing-out rotor for round bottom centrifuge tube.

SPECIFICATIONS

Rotor type	Swing-out
Tube dimensions (Ø × length)	16 × 105 mm
Capacity	12
Tube's volume	10–15 ml
Max. speed	4200 rpm
Max. RCF	3160 × g
Tube manufacturers	Nunc, Greiner, Greiner Bio-one, TPP etc.
Autoclavable	+



CAT. NUMBER

BS-010221-MK

Rotor Ri-12/10



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BN-13/75
BS-010208-PK
adapter set

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-13/75 for vacutainers 2-5 ml [13x75]



BN-13/100
BS-010208-QK
adapter set

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-13/100 for vacutainers 4-8 ml [13x100]



BN-16/100
BS-010208-RK
adapter set

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-16/100 for vacutainers 8-10 ml [16x100]



RR-U
BS-010208-UK
rack for rotors

Rack for rotors



BN-16/90
BS-010208-SK
adapter set for rotor R-12/10 R-24/10

Spare adapter set for rotor R-12/10, R-24/10. Plastic buckets (12 pcs).