

Ri-12/15, Rotor

DESCRIPTION

Swing-out rotor for conical bottom centrifuge tube.

SPECIFICATIONS

Rotor type	Swing-out
Tube dimensions (Ø × length)	17 × 120 mm
Capacity	12
Tube's volume	15 ml
Max. speed	4200 rpm
Max. RCF	3160 × g
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.
Autoclavable	+



CAT. NUMBER

BS-010221-KK

Rotor Ri-12/15

ACCESSORIES



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



RR-U
BS-010208-UK
rack for rotors

Rack for rotors



BN-17/120
BS-010208-TK
adapter set for rotor R-12/15

Spare adapter set for rotor R-12/15. Plastic buckets (12 pcs).