

Ri-2, Rotor

DESCRIPTION

Swing-out rotor for standard 96-well microtitre plates, skirted PCR plates, deepwell plates up to 45mm.

SPECIFICATIONS

Rotor type	Swing-out
Platform dimensions (w × l)	128 × 85.6 mm
Max. height	up to 45 mm
Capacity	2
Max. speed	2000 rpm
Max. RCF	560 × g
Tube manufacturers	Nunc, Greiner, Sarstedt, Corning, Greiner Bio-one, etc.
Autoclavable	+



CAT. NUMBER

BS-010221-GK

Rotor Ri-2

ACCESSORIES



LMC-56
 Without rotor
 Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



AP-96
 BS-010219-DK
 adapter

2 adapters for 96-well semi-skirted and non-skirted PCR plates - made of Ertacetal® C. Autoclavable



AP-384
 BS-010219-EK
 adapter

2 adapters for 384-well plates - made of Ertacetal® C. Autoclavable



RR-U
 BS-010208-UK
 rack for rotors

Rack for rotors