

Ri-24GC, Rotor

DESCRIPTION

Rotor for 24 Gel Cards for blood group serology testing (Forward Group, Reverse Group, RhD Type and 3 cell antibody screen)

SPECIFICATIONS

Rotor type	Swing-out
Card dimensions (l × w)	53 × 74 mm
Capacity	24 gel cards
Max. speed	1500 rpm
Max. RCF	280 × g
Gel card Manufacturers	Grifols, DiaMed, etc.
Autoclavable	+



CAT. NUMBER

BS-010221-NK

Rotor Ri-24GC

ACCESSORIES



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



RR-U
BS-010208-UK
rack for rotors

Rack for rotors