

BA-2/50, Bucket adapters

DESCRIPTION

Bucket adapters for 2 tubes of 50 ml (set of 4 pcs.)

SPECIFICATIONS

Tube holders material	Aluminium
Rotor type	Swing-out
Tube dimensions (Ø × length)	29 × 115 mm
Capacity	2
Tube's volume	50 ml
Max. speed	4200 rpm
Max. RCF	3260 × g

Tube manufacturers Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.



CAT. NUMBER

BS-010221-CK

BA-2/50

ACCESSORIES



LMC-56
 Without rotor
 Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BR-4U
 BS-010221-AK
 Rotor

Bucket rotor without bucket adapters



BI-25-8
 BS-010221-UK
 adapter set

For bucket adapters BA-2/50 optional adapter set (8 pcs) - BI-25-8 for 25 ml conical tubes



RR-U
 BS-010208-UK
 rack for rotors

Rack for rotors