

BA-4/RT, Bucket adapters

DESCRIPTION

Bucket adapters for 4 round bottom centrifuge tubes of 10 ml (set of 4 pcs.)

SPECIFICATIONS

Tube holders material	POM Kocetal
Rotor type	Swing-out
Capacity	4
Tube's volume	10–15 ml
Max. speed	4200 rpm
Max. RCF	2960 × g
Tube manufacturers	Falcon, Greiner Bio-one, Sarstedt, Corning, Nunc, etc.



CAT. NUMBER

BS-010221-3K

BA-4/RT

ACCESSORIES



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BR-4U
BS-010221-AK
Rotor

Bucket rotor without bucket adapters



RR-U
BS-010208-UK
rack for rotors

Rack for rotors