

BA-14/2U, Bucket adapters

DESCRIPTION

Bucket adapters for 14 microtest tubes of 2 ml (set of 4 pcs.)

SPECIFICATIONS

Tube holders material	POM Kocetal
Rotor type	Swing-out
Capacity	14
Tube's volume	2 ml
Max. speed	4200 rpm
Max. RCF	2070 × g
Tube manufacturers	Nunc, Greiner, Sarstedt, Corning, Greiner Bio-one, etc.



CAT. NUMBER

BS-010221-JK

BA-14/2U

ACCESSORIES



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BR-4U
BS-010221-AK
Rotor

Bucket rotor without bucket adapters



RR-U
BS-010208-UK
rack for rotors

Rack for rotors