

BR-4U, Rotor

DESCRIPTION

Bucket rotor without bucket adapters

SPECIFICATIONS

Rotor type

Swing-out

Max. speed

4200 rpm



CAT. NUMBER

BS-010221-AK

BR-4U, Rotor

ACCESSORIES



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BA-2/50
BS-010221-CK
Bucket adapters

Bucket adapters for 2 tubes of 50 ml (set of 4 pcs.)



BI-25-8
BS-010221-UK
adapter set

For bucket adapters BA-2/50 optional adapter set (8 pcs) - BI-25-8 for 25 ml conical tubes



BA-4/FT
BS-010221-YK
Bucket adapters

Bucket adapters for 4 conical bottom centrifuge tubes of 15 ml (set of 4 pcs.)



BA-4/VT
BS-010221-4K
Bucket adapters

Bucket adapters for 4 vacutainer adapters of 2-9 ml (set of 4 pcs.). Included 8-9 ml [16 x 100] vacutainer adapters, BN-16/100-16.



BN-13/75-16
BS-010221-OK
adapter set

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-13/75-16 for vacutainers 2-5 ml [13x75]



BN-13/100-16
BS-010221-PK
adapter set

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-13/100-16 for vacutainers 4-8 ml [13x100]



BN-16/100-16
BS-010221-QK
adapter set

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-16/100-16 for vacutainers 8-10 ml [16x100]



BA-4/RT
BS-010221-3K
Bucket adapters

Bucket adapters for 4 round bottom centrifuge tubes of 10 ml (set of 4 pcs.)



BA-8/5
BS-010221-FK
Bucket adapters

Bucket adapters for 8 tubes of 5 ml (set of 4 pcs.)



BA-14/2U
BS-010221-JK
Bucket adapters

Bucket adapters for 14 microtest tubes of 2 ml (set of 4 pcs.)



RR-U
BS-010208-UK
rack for rotors

Rack for rotors