

BN-13/75-16, adapter set

DESCRIPTION

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-13/75-16 for vacutainers 2-5 ml [13x75]



CAT. NUMBER

BS-010221-OK

For bucket adapters BA-4/VT optional adapter set
(16 pcs) - BN-13/75-16 for vacutainers 2-5 ml

ACCESSORIES



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BR-4U
BS-010221-AK
Rotor

Bucket rotor without bucket adapters



BA-4/VT
BS-010221-4K
Bucket adapters

Bucket adapters for 4 vacutainer adapters of 2-9 ml (set of 4 pcs.). Included 8-9 ml [16 x 100] vacutainer adapters, BN-16/100-16.