

## BN-13/100-16, adapter set

### DESCRIPTION

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-13/100-16 for vacutainers 4-8 ml [13x100]



### CAT. NUMBER

**BS-010221-PK**

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-13/100-16 for vacutainers 4-8 ml

### ACCESSORIES



**LMC-56**  
Without rotor  
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



**BR-4U**  
BS-010221-AK  
Rotor

Bucket rotor without bucket adapters



**BA-4/VT**  
BS-010221-4K  
Bucket adapters

Bucket adapters for 4 vacutainer adapters of 2-9 ml (set of 4 pcs.). Included 8-9 ml [16 x 100] vacutainer adapters, BN-16/100-16.