

BN-13/100-24, adapter set

DESCRIPTION

For rotors R-24/10, Ri-24/10 optional adapter set (24 pcs) - BN-13/100-24 for vacutainers 4-8 ml [13x100]



CAT. NUMBER

BS-010221-SK

For rotors R-24/10, Ri-24/10 optional adapter set
(24 pcs) - BN-13/100-24 for vacutainers 4-8 ml

ACCESSORIES



LMC-4200R

Without rotor
Refrigerated Laboratory
Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



LMC-56

Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



R-24/10

BS-010212-JK
Rotor

Swing-out rotor for round bottom centrifuge tube.



Ri-24/10

BS-010221-LK
Rotor

Swing-out rotor for round bottom centrifuge tube.