

BI-25-8, adapter set

DESCRIPTION

For bucket adapters BA-2/50 optional adapter set (8 pcs) - BI-25-8 for 25 ml conical tubes



CAT. NUMBER

BS-010221-UK

For bucket adapters BA-2/50 optional adapter set (8 pcs) - BI-25-8 for 25 ml conical tubes

ACCESSORIES



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BR-4U
BS-010221-AK
Rotor

Bucket rotor without bucket adapters



BA-2/50
BS-010221-CK
Bucket adapters

Bucket adapters for 2 tubes of 50 ml (set of 4 pcs.)