


B-05E, block

DESCRIPTION

Block for one deep-well plate Eppendorf® 96/500 µl - Cat.No. 0030501101



CAT. NUMBER

BS-010159-QK

B-05E block for one deep-well plate
Eppendorf® 96/500 µl



TS-DW
Without thermoblock
Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)