


B-1R, block

DESCRIPTION

Block for one deep-well plate Riplate® 96/1000 µl - Cat.No. 43001-0101



CAT. NUMBER

BS-010159-UK

Block for one deep-well plate Riplate® 96/1000
µl - Cat.No. 43001-0101



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)