


B-2R, block

DESCRIPTION

Block for one deep-well plate Riplate® 96/2000 µl - Cat.No. 43001-0103



CAT. NUMBER

BS-010159-TK

Block for one deep-well plate Riplate® 96/2000
µl



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)