


B-2SL, block

DESCRIPTION

Block for one deep-well plate Starlab® 96/1200 µl - Cat.No. E2896-0120



CAT. NUMBER

BS-010159-IK Block for one deep-well plate Starlab® 96/1200 µl



TS-DW
Without thermoblock
Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)