


B-2SD, block

DESCRIPTION

Block for one deep-well plate Slicprep™ 96 - Cat.No. AB0932



CAT. NUMBER

BS-010159-PK Block for one deep-well plate Slicprep™ 96



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)