


# B-2KF, block

## DESCRIPTION

Block for one deep-well plate PowerMag® Glass Bead Plate 96 - Cat.No. 27600-4-KF-BP



## CAT. NUMBER

**BS-010159-LK**

Block for one deep-well plate PowerMag® Glass  
Bead Plate 96



### TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)