

LTD-1, Drawer

DESCRIPTION

Drawer with non slip silicone mat and soft close function. Can be mounted afterwards.

SPECIFICATIONS

Drawer internal dimensions (W×D×H)	310 × 350 × 40 mm
Overall dimensions (W×D×H)	430 × 410 × 70 mm



CAT. NUMBER

BS-040107-OK

LTD-1, Drawer

ACCESSORIES



LT-120
BS-040107-EK
Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)



LT-150
BS-040107-JK
Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)



LT-180
BS-040107-HK
Height adjustable desk

This laboratory table is the perfect solution for any research or laboratory facility. It features a chemically resistant high pressure laminate worktop, making it suitable for use in a variety ...

[read more](#)