

# Storage tank 60l

## DESCRIPTION

Storage tank with base, tap and multipoint level switch, 60 l



## CAT. NUMBER

BS-070102-SK

Storage tank, 60l

## ACCESSORIES



**Labaqua Trace**  
BS-070105-A02  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labaqua ultrapure system produces both ultrapure ...

[read more](#)



**Labaqua HPLC**  
BS-070104-A02  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

**Labaqua HPLC** produces water with very low organic carbon (TOC) ...

[read more](#)



**Labaqua Bio**  
30l storage tank  
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

**Labaqua Bio** system produces water with very low organic and ...

[read more](#)