

CO₂ Incubator

DESCRIPTION

The CO₂ incubator is the most important piece of equipment for biosafety laboratory R&D. The Haier Biomedical CO₂ incubator uses professional and superior technology to provide reliable conditions for frontline medical staff to carry out epidemiological and translational research, and to create a stable cell and virus culture environment for researchers.

- Efficient microbial filter
- Active airflow control technology, to ensure stability and homogeneity
- Door heating function, eliminates microbial contamination
- 180°C dry heat sterilization technology, providing a contamination-free environment

[More info](#)

